This package includes (1) Foot File, (1) protective cover, and (1) skin receptacle.

The Microplane® Foot File:

- is intended for removing dry skin and calluses on the heels and soles of the feet
- is specifically designed to work BEST on DRY skin, no need for pre-soaking
- can be sanitized or steam sterilized following recommended procedures – skin receptacle and protective cover may only be sanitized in Barbicide® solution

Foot File Use:

1. For best results, use on dry, moisture-free skin.
2. Install skin receptacle into backside of file surface with smaller end located nearest to the handle.
3. Remove protective cover from face of file.
4. Gently move Foot File in a back and forth motion (in direction of rounded file surface ends) across the area of the foot to be treated, adding slightly more pressure as needed to achieve the desired result. AVOID moving the Foot File in a side-to-side or circular motion.
5. Stop use when new skin is revealed or when callus has been removed, as overuse may cause feet to become sore.
6. To dispose of skin filings, remove skin receptacle from Foot File by pressing on lined press point with thumb. The opposite end will pop up slightly: Use fingertip to loosen.
7. Clean Foot File and skin receptacle after each use (see below).
8. Secure protective cover over file surface and re-install skin receptacle for storage between uses.

Cleaning Procedures:

Sanitizing Process – Foot File

1. Remove skin receptacle and protective cover from Foot File.
2. Flush front of file surface with running water for 5 seconds while slowly tilting paddle from end to end (seesaw action).
3. Turn file over and flush back of file surface with running water for 5 seconds, repeating the tilting motion described in step 2.
4. Prepare proper mixture of Barbicide® solution by following directions on Barbicide® container for “Beauty/Barber Instruments and Tools.”
5. Dip head of file completely in solution, insuring that the entire metal file surface is totally submerged. Move file head in a back-and-forth motion in the solution for 5 seconds.
6. Submerge the entire Foot File in the solution and allow it to soak for 10 minutes.
7. Remove Foot File from solution. Using tilting motion described in step 2, flush both the front and back of file surface with running water for 5 seconds each in order to remove any skin particles that might have settled on the tool while soaking.
8. Allow to air dry before next use.

Sanitizing Process – Protective Cover and Skin Receptacle

1. Flush both the protective cover and skin receptacle with running water for 5 seconds (front and back).
2. Prepare proper mixture of Barbicide® solution by following directions on Barbicide® container for “Beauty/Barber Instruments and Tools.”
3. Submerge pieces in solution and allow to soak for 10 minutes.
4. Remove from solution and flush them with running water for 5 seconds (front and back).

5. Allow to air dry before next use.

Sterilization Process – Foot File only

1. Remove skin receptacle and protective cover from Foot File.
2. Flush front of file surface with running water for 5 seconds while slowly tilting paddle from end to end (seesaw action).
3. Turn file over and flush back of file surface with running water for 5 seconds, repeating the tilting motion described in step 2.
4. Tightly wrap entire Foot File in aluminum foil and place in steam autoclave at 275 °F, 15 psi for 15 minutes.
5. Following autoclave manufacturer’s instructions, precautions, and recommended personal protective equipment, remove sterilized parts from autoclave.
6. Allow to air dry before next use.

CAUTION: Do not use this product if you are diabetic, have poor blood circulation, or have any inflammation or open wounds on your feet. Consult your physician before using.

Microplane®...from the woodshop...to the kitchen...to the salon...an evolution to foot care revolution

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